

Product Datasheet

anti- TriMethyl-Histone H3-K27 antibody (orb1784841)

Description	anti- TriMethyl-Histone H3-K27 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic trimethylated peptide around K27 of human histone H3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Purity	≥95% as determined by SDS-PAGE
MW	16 kDa
Uniprot ID	Q16695
Dilution Range	WB : 1:500-1:1000
Expiration Date	12 months from date of receipt.