

Product Datasheet

18:1 TAP (DOTAP, MS Salt) (orb1785909)

Description	18:1 TAP is a cationic liposome-forming compound that may be used for transfection of DNA, RNA, and ...
Note	For research use only
Formula	$C_{42}H_{80}NO_4 \cdot CH_3O_4S^-$
MW	774.18602
CAS Number	[144189-73-1]
SMILES	<chem>C(C(OC(CCCCCC/C=C/CCCCCCC)=O)COC(CCCCCC/C=C/CCCCCCC)=O)[N+](C)(C)C.COS([O-])(=O)=O</chem>
Expiration Date	12 months from date of receipt.