


Product Datasheet

Human PIGR ELISA Kit (orb1786040)

Description	Human PIGR ELISA Kit (96 Tests). Quantitate Human PIGR in cell culture supernatants, serum and plasm...
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K19-R638
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	P01833
Sensitivity	150 pg/ml
Expiration Date	Please enquire.



Human PIGR
PicoKine
ELISA Kit
standard c...