


Product Datasheet

Anti-Alpha 1 microglobulin/Ambp Antibody (orb1786069)

Description	Anti-Alpha 1 microglobulin/Ambp Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	E.coli-derived mouse Alpha 1 microglobulin/Ambp recombinant protein (Position: A205-S349).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Western blot, 0.25-0.5 µg/ml, Mouse, Rat ELISA, 0.1-0.5 µg/ml, Mouse. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q07456
Expiration Date	12 months from date of receipt.



WB analysis of
Alpha 1
Microglobulin/Ambp...