



Product Datasheet

Anti-NUBP1 Antibody (orb1786114)


Description	Anti-NUBP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	E.coli-derived human NUBP1 recombinant protein (Position: M1-Q311).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Western blot, 0.25-0.5 µg/ml, Human Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry (Fixed), 1-3 µg/1x10 ⁶ cells, Human ELISA, 0.1-0.5 µg/ml, Human. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35 kDa
Uniprot ID	P53384
Expiration Date	12 months from date of receipt.



Flow
Cytometry
analysis of
K562 cells
us...



WB analysis of
NUBP1 using
anti-NUBP1
an...



IF analysis of
NUBP1 using
anti-NUBP1
an...