

Product Datasheet

**Mouse NKp46 / NCR1 / CD335 Protein, Fc Tag
(orb1786128)**

Description

Mouse NKp46 Protein, Fc Tag is

(HEK293). It contains AA Glu 22 - Asn
...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

NKp46

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 52.9 kDa.

Purity

90%

Source

Mouse NKp46 Protein, Fc Tag (orb1786128) is expressed from human 293 cells (HEK293). It contains AA Glu 22 - Asn 255 (Accession # Q8C567-1).

MW

52.9 kDa

Uniprot ID

Q8C567

Expiration Date

6 months from date of receipt.