

Product Datasheet

Human LPA Protein, His Tag (orb1786131)

Description

Human LPA Protein, His Tag is expressed

contains AA Glu 20 - Asn 2...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LPA

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 226.3 kDa. The protein migrates as > 250 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human LPA Protein, His Tag (orb1786131) is expressed from human 293 cells (HEK293). It contains AA Glu 20 - Asn 2040 (Accession # P08519).

MW

226.3 kDa

Uniprot ID
P08519
Expiration Date

6 months from date of receipt.