
Product Datasheet


**Human IGFBP-1 Protein, His Tag (MALS verified)
(orb1786132)**

Description

Human IGFBP-1 Protein, His Tag is expressed

Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IGFBP-1
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 27.2 kDa. The protein migrates as 32-35 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human IGFBP-1 Protein, His Tag (orb1786132) is expressed from human 293 cells (HEK293). It contains AA Ala 26 - Asn 259 (Accession # P08833).
MW	27.2 kDa
Uniprot ID	P08833
Expiration Date	6 months from date of receipt.

Human IGFBP-1 Protein, His Tag (100)


 Immobilized
Human IGF-
I, Fc Tag
(Cat. No...)