

## Product Datasheet

### Mouse IL-21 R Protein, Fc Tag (orb1786142)

**Description** Mouse IL-21 R Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Cys 20 - Pr...

**Reactivity** Mouse

**Endotoxins** 1.0 EU per µg

**Conjugation** Unconjugated

**Target** IL-21 R

**Preservatives** Powder

**Form/Appearance** Powder

**Storage** -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

**Tag** C-Fc

**Note** For research use only

**Application notes** This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.5 kDa.


**Purity** 90%

**Source** Mouse IL-21 R Protein, Fc Tag (orb1786142) is expressed from human 293 cells (HEK293). It contains AA Cys 20 - Pro 236 (Accession # Q9JHX3).

**MW** 51.5 kDa

**Uniprot ID** [Q9JHX3](#)

**Expiration Date** 6 months from date of receipt.



Immobilized  
Mouse IL-21  
Protein, His  
Tag...