

Product Datasheet

Human RPE65 Protein, His Tag (orb1786145)

Description	Human RPE65 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 2 - Ser ...
Reactivity	Human
Conjugation	Unconjugated
Target	RPE65
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 62.8 kDa.
Source	Human RPE65 Protein, His Tag (orb1786145) is expressed from human 293 cells (HEK293). It contains AA Ser 2 - Ser 533 (Accession # Q16518).
MW	62.8 kDa
Uniprot ID	Q16518
Expiration Date	6 months from date of receipt.