

---

## Product Datasheet

**Human Follistatin Protein, His Tag, premium grade  
(orb1786154)**

**Description**

Human Follistatin Protein, His Tag,

human 293 cells (HEK293). It con...

**Reactivity**

Human

**Endotoxins**

0.1 EU per µg

**Conjugation**

Unconjugated

**Target**

Follistatin

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 33.4 kDa . The protein migrates as 35-36 kDa, 38-40 kDa and 41-43 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Human Follistatin Protein, His Tag, premium grade (orb1786154) is expressed from human 293 cells (HEK293). It contains AA Gly 30 - Asn 317 (Accession # P19883-2).

**MW**

33.4 kDa

**Uniprot ID****P19883****Expiration Date**

6 months from date of receipt.