

Product Datasheet

**Human TrkC / NTRK3 Protein, Fc Tag (MALS
verified) (orb1786159)**

Description

Human TrkC Protein, Fc Tag is

(HEK293). It contains AA Cys 32 - Thr 4...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TrkC

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 71.2 kDa. The protein migrates as 95-130 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human TrkC Protein, Fc Tag (orb1786159) is expressed from human 293 cells (HEK293). It contains AA Cys 32 - Thr 429 (Accession # Q16288-1).

MW

71.2 kDa

Expiration Date

6 months from date of receipt.