

Product Datasheet

Human PCSK9 Protein, His Tag (orb1786162)

Description

Human PCSK9 Protein, His Tag is

(HEK293). It contains AA Gln 31 -His ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

PCSK9

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 13.8 kDa & 33.5 kDa. The protein migrates as 15 kDa and 33 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human PCSK9 Protein, His Tag (orb1786162) is expressed from human 293 cells (HEK293). It contains AA Gln 31 -His 449 (Accession # Q8NBP7-1).

MW

13.8 kDa & 33.5 kDa

Uniprot ID

Q8NBP7

Expiration Date

6 months from date of receipt.