

Product Datasheet

Human CEACAM-5 / CD66e (499-685) Protein, His Tag (MALS verified) (orb1786166)

Description

Human CEACAM-5 (499-685)

from human 293 cells (HEK293).
It contains AA...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CEACAM-5

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in
lyophilized state;. -70°C for 3
months under sterile conditions
after reconstitution. For long
term storage, the product
should be stored at lyophilized
state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a
polyhistidine tag at the C-
terminus. The protein has a
calculated MW of 21.9 kDa. The
protein migrates as 40-55 kDa
when calibrated against Star
Ribbon Pre-stained Protein
Marker under reducing (R)
condition (SDS-PAGE) due to
glycosylation.

Purity

● -20°C to -70°C for 12 months
in lyophilized state; ● -70°C for
3 months under sterile
conditions after reconstitution.
For long term storage, the
product should be stored at
lyophilized state at -20°C or
lower.

Source

Human CEACAM-5 (499-685)
Protein, His Tag (orb1786166) is
expressed from human 293
cells (HEK293). It contains AA
Glu 499 - Ala 685 (Accession #
P06731-1).