

Product Datasheet

**Mouse LIV-1 / SLC39A6 Protein, Fc Tag (MALS
verified) (orb1786169)**

Description

Mouse LIV-1 Protein, Fc Tag is

(HEK293). It contains AA Leu 21 - Trp ...

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LIV-1

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 62.3 kDa. The protein migrates as 80-95 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse LIV-1 Protein, Fc Tag (orb1786169) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Trp 335 (Accession # Q8C145-1).

MW

62.3 kDa

Uniprot ID

Q8C145

Expiration Date

6 months from date of receipt.