

Product Datasheet

**Human OX40 Ligand / TNFSF4 Protein, His Tag,
active trimer (MALS verified) (orb1786170)**

Description

Human OX40 Ligand Protein, His

293 cells (HEK293). It contains AA Gln 51...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

OX40 Ligand

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 49.4 kDa. The protein migrates as 50-116 kDa under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.

Purity

80%

Source

Human OX40 Ligand Protein, His Tag (orb1786170) is expressed from human 293 cells (HEK293). It contains AA Gln 51 - Leu 183 (Accession # P23510-1).

MW

49.4 kDa

Uniprot ID
P23510
Expiration Date

6 months from date of receipt.