

## Product Datasheet

**Rhesus macaque TROP-2 / TACSTD2 Protein, Fc Tag (MALS verified) (orb1786171)**

**Description**

Rhesus macaque TROP-2 Protein,

293 cells (HEK293). It contains  
AA Gln...

**Reactivity**

Monkey

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

TROP-2

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in  
lyophilized state;. -70°C for 3  
months under sterile conditions  
after reconstitution. For long  
term storage, the product should  
be stored at lyophilized state at -  
20°C or lower

**Tag**

C-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1  
Fc tag at the C-terminus. The  
protein has a calculated MW of  
53.9 kDa. The protein migrates  
as 65-80 kDa when calibrated  
against Star Ribbon Pre-stained  
Protein Marker under reducing  
(R) condition (SDS-PAGE) due to  
glycosylation.

**Purity**

90%

**Source**

Rhesus macaque TROP-2 Protein,  
Fc Tag (orb1786171) is  
expressed from human 293 cells  
(HEK293). It contains AA Gln 31 -  
Thr 274 (Accession #  
XP\_001114599.1).

**MW**

53.9 kDa

**Uniprot ID**

XP\_001114599.1

**Expiration Date**

6 months from date of receipt.