

Product Datasheet

Human TROP-2 / TACSTD2 (88-274) Protein, His Tag (orb1786173)

Description

Human TROP-2 (88-274) Protein,

293 cells (HEK293). It contains AA Thr...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TROP-2

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.5 kDa. The protein migrates as 33-43 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human TROP-2 (88-274) Protein, His Tag (orb1786173) is expressed from human 293 cells (HEK293). It contains AA Thr 88 - Thr 274 (Accession # P09758-1).

MW

23.5 kDa

Uniprot ID
P09758
Expiration Date

6 months from date of receipt.