

## Product Datasheet

**Human Transferrin R / CD71 Protein, Fc Tag (MALS  
verified) (orb1786176)**

**Description**

Human Transferrin R Protein, Fc

cells (HEK293). It contains AA Cys 8...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

Transferrin R

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 101.6 kDa.

**Purity**

95%

**Source**

Human Transferrin R Protein, Fc Tag (orb1786176) is expressed from human 293 cells (HEK293). It contains AA Cys 89 - Phe 760 (Accession # P02786-1).

**MW**

101.6 kDa

**Uniprot ID**
**P02786**
**Expiration Date**

6 months from date of receipt.