

Product Datasheet

**Human TGF-beta RI Protein, Fc Tag (MALS
verified) (orb1786177)**

Description

Human TGF-beta RI Protein, Fc Tag

(HEK293). It contains AA Leu 34 ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TGF-beta RI

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 36.7 kDa.

Purity

90%

Source

Human TGF-beta RI Protein, Fc Tag (orb1786177) is expressed from human 293 cells (HEK293). It contains AA Leu 34 - Leu 126 (Accession # P36897-1).

MW

36.7 kDa

Uniprot ID
P36897
Expiration Date

6 months from date of receipt.