

Product Datasheet

Cynomolgus CD6 Protein, His Tag (MALS verified) (orb1786181)

Description

Cynomolgus CD6 Protein, His Tag is

(HEK293). It contains AA His 18 - ...

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CD6

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 42.8 kDa. The protein migrates as 55-80 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Cynomolgus CD6 Protein, His Tag (orb1786181) is expressed from human 293 cells (HEK293). It contains AA His 18 - Leu 399 (Accession # XP_005577746.1).

MW

42.8 kDa

Uniprot ID

XP_005577746.1

Expiration Date

6 months from date of receipt.