

Product Datasheet

Rabies virus (strain ERA) Nucleoprotein, Tag Free (orb1786190)

Description

Rabies virus (strain ERA)

expressed from human 293 cells (HEK293). It con...

Reactivity

Rabies virus

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Nucleoprotein

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 50.7 kDa. The protein migrates as Band dispersion when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

70%

Source

Rabies virus (strain ERA)
Nucleoprotein, Tag Free (orb1786190) is expressed from human 293 cells (HEK293). It contains AA Met 1 - Ser 450 (Accession # P0DOF3).

MW

50.7 kDa

Uniprot ID
P0DOF3
Expiration Date

6 months from date of receipt.