

---

## Product Datasheet

**Human CCL20 Protein, His Tag (HPLC verified)  
(orb1786198)**

**Description**

Human CCL20 Protein, His Tag is

contains AA Ala 27 - Met 96  
(Access...**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

CCL20

**Preservatives**

0.085% TFA 30% ACN

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 10 kDa. The protein migrates as 11-12 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

**Purity**

90%

**Source**

Human CCL20 Protein, His Tag (orb1786198) is expressed from E. coli cells. It contains AA Ala 27 - Met 96 (Accession # P78556-1).

**MW**

10 kDa

**Uniprot ID****P78556****Expiration Date**

6 months from date of receipt.