
Product Datasheet

**Human Progranulin Protein, His Tag, low
endotoxin (MALS & SPR verified) (orb1786203)**

Description

Human Progranulin Protein, His

293 cells (HEK293). It contains AA Thr 18...

Reactivity

Human

Endotoxins

0.01 EU per µg

Conjugation

Unconjugated

Target

Progranulin

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 63.7 kDa. The protein migrates as 80-95 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human Progranulin Protein, His Tag (orb1786203) is expressed from human 293 cells (HEK293). It contains AA Thr 18 - Leu 593 (Accession # P28799-1).

MW

63.7 kDa

Uniprot ID**P28799****Expiration Date**

6 months from date of receipt.