
Product Datasheet

**Human FGF basic (154aa) Protein, premium grade
(orb1786204)**

Description

Human FGF basic (154aa) Protein,

coli cells. It contains AA Ala 1...

Reactivity

Human

Endotoxins

0.01 EU per µg

Conjugation

Unconjugated

Target

FGF basic

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 17.3 kDa. The protein migrates as 15-17 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

Purity

95%

Source

Human FGF basic (154aa) Protein, premium grade (orb1786204) is expressed from E. coli cells. It contains AA Ala 135 - Ser 288 (Accession # NP_001997.5).

MW

17.3 kDa

Uniprot ID

NP_001997

Expiration Date

6 months from date of receipt.