

Product Datasheet

Anti-HA-tag Antibody-HRP labeled (orb1786243)

Description	Rabbit polyclonal antibody to HA-tag-HRP labeled
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of HA-tag.
Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Proclin
Concentration	1 mg/ml
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Dilution Range	WB (1/2000 - 1/5000), IC (1/500 - 1/2000)
Expiration Date	12 months from date of receipt.