

## Product Datasheet

### Goat Anti-Rabbit IgG (H&L) - Cy7 (orb1786294)

<b>Description</b>	Goat Polyclonal Secondary Antibody to Rabbit IgG (H&L) Cy7 labeled
<b>Species/Host</b>	Goat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC, IF, IH
<b>Immunogen</b>	Rabbit IgG
<b>Preservatives</b>	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Dilution Range</b>	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
<b>Expiration Date</b>	12 months from date of receipt.