

Product Datasheet

Rabbit Anti-Goat IgG (H&L) - AP (orb1786498)

| | |
|----------------------------|--|
| Description | Rabbit Polyclonal Secondary Antibody to Goat IgG (H&L) AP labeled |
| Species/Host | Rabbit |
| Conjugation | AP |
| Tested Applications | ELISA, IH, WB |
| Immunogen | Goat IgG |
| Preservatives | 0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide |
| Concentration | 1 mg/ml |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Dilution Range | E (1/10000 - 1/50000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000) |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)