

Product Datasheet

Goat Anti-Human IgG F(ab')₂ - Texas Red (orb1786525)

Description	Goat Polyclonal Secondary Antibody to Human IgG F(ab') ₂ Texas Red labeled
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IH
Immunogen	Human IgG F(ab') ₂
Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Concentration	1 mg/ml
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Dilution Range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
Expiration Date	12 months from date of receipt.