

## Product Datasheet

### Goat Anti-Mouse IgM - TRITC (orb1786574)

|                            |                                                                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat Polyclonal Secondary Antibody to Mouse IgM TRITC labeled                                                                                    |
| <b>Species/Host</b>        | Goat                                                                                                                                             |
| <b>Conjugation</b>         | TRITC                                                                                                                                            |
| <b>Tested Applications</b> | FC, IF, IH                                                                                                                                       |
| <b>Immunogen</b>           | Mouse IgM                                                                                                                                        |
| <b>Preservatives</b>       | 0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                                          |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                                 |
| <b>Note</b>                | For research use only                                                                                                                            |
| <b>Dilution Range</b>      | IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000)                                                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                  |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)