

## Product Datasheet

# Rabbit Anti-Donkey IgG (H&L) - TRITC (orb1786602)

<b>Description</b>	Rabbit Polyclonal Secondary Antibody to Donkey IgG (H&L) TRITC labeled
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	TRITC
<b>Tested Applications</b>	ELISA, FC, IF, IH
<b>Immunogen</b>	Donkey IgG
<b>Preservatives</b>	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Dilution Range</b>	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)