

Product Datasheet

Rabbit Anti-Donkey IgG (H&L) - FITC (orb1786705)

Description	Rabbit Polyclonal Secondary Antibody to Donkey IgG (H&L) FITC
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Donkey IgG
Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Concentration	1 mg/ml
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Dilution Range	IF (1/250 - 1/2000), FC (1/250 - 1/2000)
Expiration Date	12 months from date of receipt.