

Product Datasheet

Rabbit Anti-Goat IgG (H&L) - Biotin (orb1786759)

Description	Rabbit Polyclonal Secondary Antibody to Goat IgG (H&L) Biotin
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IH, WB
Immunogen	Goat IgG
Preservatives	0.02% Proclin/Sodium Azide. 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Sodium Azide
Concentration	1 mg/ml
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Dilution Range	E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.