


Product Datasheet

Mouse MIP-3 Alpha/CCL20 ELISA Kit (orb178676)

Description	Mouse MIP-3 Alpha/CCL20 ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A28-M97
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 0-26 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10826 MW
Uniprot ID	089093
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse MIP-3
alpha/CCL20
PicoKine
ELISA K...