


## Product Datasheet

### Human CXCL7/PPBP/NAP-2 ELISA Kit (orb178686)

<b>Description</b>	Human CXCL7/PPBP/NAP-2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6 pg/ml - 1,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: A59-D128
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:500, concentration in serum and plasma is 50-150 ng/ml
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	13894 MW
<b>Uniprot ID</b>	<b>P02775</b>
<b>Sensitivity</b>	2 pg/ml
<b>Expiration Date</b>	Please enquire.



Human  
CXCL7  
PicoKine  
ELISA Kit  
standard ...