


Product Datasheet

Human Kallikrein-6 / KLK6 ELISA Kit (orb178688)

Description	Human Kallikrein-6 / KLK6 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L22-K244
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	26856 MW
Uniprot ID	Q92876
Sensitivity	3 pg/ml
Expiration Date	Please enquire.



Human Kallikrein-6 PicoKine ELISA Kit st...