


Product Datasheet

Human NOV/CCN3 ELISA Kit (orb178689)

| | |
|----------------------------|---|
| Description | Human NOV/CCN3 ELISA Kit |
| Reactivity | Human |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: T32-M357 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| Clonality | Polyclonal |
| MW | 39162 MW |
| Uniprot ID | P48745 |
| Sensitivity | 2 pg/ml |
| Expiration Date | Please enquire. |



Human NOV/CCN3 PicoKine ELISA Kit standa...