

Product Datasheet

JNK1/2/3 (Phospho-T183/Y185) polyclonal antibody \square KU \square (orb1786934)

Description	JNK1/2/3 (Phospho-T183/Y185) polyclonal antibody \square KU \square
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P53779 , P45983 , P45984
Expiration Date	12 months from date of receipt.