


Product Datasheet

Mouse CCL9/MIP-1 gamma ELISA Kit (orb178700)

Description	Mouse CCL9/MIP-1 gamma ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q22-Q122
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1000, concentration in serum and plasma is 50-200 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	13871 MW
Uniprot ID	P51670
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Mouse CCL9
PicoKine
ELISA Kit
standard
C...