

Product Datasheet

Symmetric DiMethyl-Histone H3-R8 polyclonal antibody (orb1787329)

Description	Symmetric DiMethyl-Histone H3-R8 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IF, IHC, IP, WB
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q16695
Expiration Date	12 months from date of receipt.