

Product Datasheet

α -Smooth Muscle Actin (ACTA2) polyclonal antibody (orb1787538)

Description	α -Smooth Muscle Actin (ACTA2) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human α -Smooth Muscle Actin (ACTA2)(NP_001135417.1).
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	42KDa
Uniprot ID	P62736
Dilution Range	WB,1:500 - 1:1000 IHC,1:50 - 1:100 IF/ICC,1:100 - 1:200
Expiration Date	12 months from date of receipt.