

Product Datasheet

p70 S6 kinase α (Phospho-Thr229) polyclonal antibody \square KU \square (orb1787579)

| | |
|----------------------------|--|
| Description | p70 S6 kinase α (Phospho-Thr229) polyclonal antibody \square KU \square |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Concentration | 1mg/ml |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | P23443 |
| Expiration Date | 12 months from date of receipt. |