

## Product Datasheet

# Smad2 (Phospho-Ser250) polyclonal antibody [KU] (orb1787580)

<b>Description</b>	Smad2 (Phospho-Ser250) polyclonal antibody [KU]
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">Q15796</a>
<b>Expiration Date</b>	12 months from date of receipt.