

Product Datasheet

Akt (Phospho-Ser124) polyclonal antibody [KU] (orb1787582)

Description	Akt (Phospho-Ser124) polyclonal antibody [KU]
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P31749 , P31751 , Q9Y243
Expiration Date	12 months from date of receipt.