

## Product Datasheet

# Recombinant His Tag Antibody / Histidine [r6HIS/6423] (orb1787777)

<b>Description</b>	Recombinant His Tag Antibody / Histidine [r6HIS/6423]
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, WB
<b>Immunogen</b>	Hexa-histidine tagged human CD70 recombinant protein was used as the immunogen for the recombinant His Tag antibody.
<b>Preservatives</b>	1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.2 mg/ml
<b>Storage</b>	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Western blot : 1-2ug/ml Flow cytometry : 1-2ug/million cells Immunofluorescence : 1-2ug/ml
<b>Isotype</b>	IgG1, kappa
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.