
Product Datasheet


StAR Antibody / Steroidogenic acute regulatory protein (orb1787814)

Description


Steroidogenic Acute Regulatory Protein (STAR)




	by ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 39-108 from the human protein was used as the immunogen for this StAR antibody.
Preservatives	1 mg/ml in 1X PBS; rAlbumin free, sodium azide free
Storage	Store the StAR antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the StAR antibody to be titered up or down for optimal performance.
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P49675
Hazard Information	This StAR antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human adrenal gland...



IHC staining of FFPE human testicular ca...



IHC staining of FFPE human testicular ca...