

Product Datasheet

Glu-Glu tag Antibody (orb1787850)

Description	Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Target	This Glu-Glu tag antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide of Glu-Glu tag.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42 Da
Uniprot ID	NA
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.