

Product Datasheet

GST Tag Antibody (orb1787853)

Description	Purified Mouse Monoclonal Antibody (Mab)
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	Purified recombinant GST fusion protein was used to produced this monoclonal antibody.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	45-50KD
Dilution Range	WB: 1:1000, WB: 1:1000
Expiration Date	12 months from date of receipt.