

## Product Datasheet

### FLAG tag antibody (DYKDDDDK) (orb1787857)

**Description** Purified Rabbit Polyclonal Antibody (Pab)

**Species/Host** Rabbit

**Conjugation** Unconjugated

**Tested Applications** WB

**Target** KLH conjugated synthetic peptide was used as antigen.

**Storage** Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

**Note** For research use only

**Isotype** Rabbit IgG

**Clonality** Polyclonal

**MW** 45-50KD

**Dilution Range** WB: 1:1000

**Expiration Date** 12 months from date of receipt.