

Product Datasheet

FEZ2 Antibody (N-term) (orb1787927)

Description	Affinity Purified Rabbit Polyclonal Antibody (Pab)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	1-30 aa
Target	This FEZ2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1-30 amino acids from the N-terminal region of human FEZ2.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39666 Da
Uniprot ID	Q9UHY8
NCBI	NP_001036013.1 , NP_005093.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.